LungPoint - Navigation and Planning System

Initiative Type
Technology
Status
Close
Added
20 September 2018
Last updated
13 June 2019
URL
https://clinicalexcellence.qld.gov.au/improvement-exchange/lungpoint-navigation-and-planning-system
Summary

Pathological confirmation of lung cancer enables definitive diagnosis of cancer and differentiates the two major types of lung cancer which are managed quite differently. The choice of diagnostic technique to diagnose lung cancer depends on the performance of each technique, the likely type and size of lung cancer and overall health of the patient. LungPoint is a new diagnostic tool

Satellite Planning System. The main goals are to maximize the yield for both diagnostic and staging purposes and to minimise risk to the patient.
Key dates
Jan 2014
Oct 2017
Implementation sites
The Prince Charles Hospital
Key Contacts
Jacqui Thomson
8160
<u>Anonymous</u>
Manager
Healthcare Improvement Unit
(07) 3328 9174
secretariat_hta@health.qld.gov.au
Aim

Provides an opportunity to pilot and evaluate new technologies within 'real world' clinical settings in

comprised of two components: 1. Virtual Bronchoscopic Navigation (VBN) System and 2. LungPoint

the Queensland context.

Benefits

Benefits of LungPoint include:

- In a small number of cases, LungPoint contributed to improved bronchoscopic navigation to the lesion by experienced Bronchoscopists.
- LungPoint allowed collaborative teaching and discussion with trainees.

Background

This technology was funded through the New Technology Funding and Evaluation Program (NTFEP). The NTFEP funds the introduction and evaluation of new technologies that:

- Are safe and effective
- Provide better health outcomes
- Provide value for money
- Provide greater access to care.

The evaluation findings will inform recommendations regarding the future use and/or investment of the technology within Queensland.

Evaluation and Results

- Beneficial to improved bronchoscopic navigation to the lesion in a small number of cases. LungPoint pathways would often result in a significantly lower diagnostic yield than EBUS RP.
- In many cases, the best path determined by the operator using 3-view CT scans was different from the path provided by Lungpoint.

Lessons Learnt

 Planning of pathways took an additional 5 minutes (approximatly) compared to current processes.

User-defined pathways are not permitted.
Resources
Technology Evaluation Summary
PDF saved 03/05/2024